

Sodium Citrate Fill Guide

BD Vacutainer® Plus Plastic Citrate Tube

BD Vacutainer® Plus Plastic Citrate Tube Draw Volume Guide

*Sufficient volume achieved if blood drawn falls above minimum fill indicator. For blood transfer, **do not** fill above illustrated dashed maximum line.*

Note: The quantity of blood drawn into evacuated tubes varies with altitude, ambient temperature, barometric pressure, tube age, venous pressure and filling technique.

Ref. # 363083

2.7 mL Draw Tube
13 mm x 75 mm Full Draw

Ref. # 363080

1.8 mL Draw Tube
13 mm x 75 mm Full Draw

Maximum Fill*

Minimum Fill Indicator

Now with Minimum Fill Indicator

Represents minimum volume of blood required for appropriate analysis

*According to CLSI guideline, Dec. 2003, Doc. H1-A5, Vol. 23, No. 33.