


SWAB DEVICES FOR MICRO AND MICRO MOLECULAR TESTS

TEST	DESCRIPTION	CONTAINER(S)	ALTERNATIVE IF AVAILABLE
LAB3513	Bacterial Vaginosis Exam	BD E-swab (White Cap)	
LAB1182	Bordetella Pertussis/Parapertussis DNA Detection	BD VTM (Red Cap, Pink Liquid)	
LAB1114	Chlamydia/GC by PCR	Alinity m multi-collect kit	
LAB2890	Covid-19 Coronavirus PCR	BD VTM (Red Cap, Pink Liquid)	
LAB1850	Covid-19 Influenza A/B RSV PCR Assay	BD VTM (Red Cap, Pink Liquid)	
LAB6699	Covid-19 Influenza A/B LIAT PCR Assay	BD VTM (Red Cap, Pink Liquid)	
LAB35	Culture, Aerobic/Anaerobic	Anaerobic Transport Medium (ATM)	BD E-swab (White Cap)
LAB34	Culture, Aerobic	BD E-swab (White Cap)	
LAB3160	Culture, Group B Strep, Prenatal Screen	BD E-swab (White Cap)	
LAB37	Culture, GC Screen	BD E-swab (White Cap)	
LAB77	Culture, Genital	BD E-swab (White Cap)	
LAB3702	Culture, Nasal	BD E-swab (White Cap)	
LAB3521	Culture, Respiratory, Cystic Fibrosis	BD E-swab (White Cap)	Sterile container
LAB2771	Culture, Rectal, Pre-Prostate Biopsy	BD E-swab (White Cap)	
LAB393	Culture, Throat Strep (only)	BD E-swab (White Cap)	
LAB3514	Culture, Yeast, Vaginal/Oral	BD E-swab (White Cap)	
LAB2845	Genital Mycoplasmas NAA, Urine	APTIMA(R) COMBO 2 urine	
LAB1123	Herpes Simplex 1 & 2 (HSV) DNA Detection	BD VTM (Red Cap, Pink Liquid)	CSF in sterile container
LAB3988	MRSA Admission Screen	BD E-swab (White Cap)	M-Swab Collection (Blue Cap)
LAB2843	Mycoplasma Hominis DNA Detection, PCR	APTIMA(R) Unisex Endo/Urethral swab	
LAB3900	Pre-Surgical Staphylococcus Screen	BD E-swab (White Cap)	M-Swab Collection (Blue Cap)
LAB6698	Respiratory Pathogen Panel with COVID-19	BD VTM (Red Cap, Pink Liquid)	
LAB1864	Vaginosis Panel NAA	BD Molecular Swab Collection Kit	
LAB1155	Varicella-Zoster Virus (VZV) DNA Detection	BD VTM (Red Cap, Pink Liquid)	CSF in sterile container

CONTAINER	SWAB DEVICE
<p>Anaerobic Transport Medium (ATM) Place specimen on top of gel (<u>do not</u> push specimen into gel)</p>	
<p>BD E-swab (White Cap)</p>	
<p>BD Molecular Swab Collection Kit</p>	
<p>BD VTM (Red Cap, Pink Liquid)</p>	
<p>APTIMA(R) Unisex – Endocervical and Male Urethral Swab White swab is a cleaning swab and not be used for collection. The blue-shafted swab labeled “unisex collection swab” used for specimen collection.</p>	
<p>M-swab Collection (Blue Cap)</p>	
<p>APTIMA(R) COMBO 2 urine collection</p>	
<p>Abbott Alinity m multi-Collect Specimen Collection Kit: Orange top tube w/swab & transfer pipette</p>	