





Department of Pathology and Laboratory Medicine

Urine Collection Device Grid




Version 08/2025

Test Name	SoftLab Code	Epic Code	Label Code	Collection Device
Urinalysis Complete Macroscopic Microscopic	UAC1 UMAC UMIC1	6620 6650 6660	UAP UAP UAP	Red/Yellow Tiger Top Tube 
24 hour Urine Write the volume of 24hr urine containers provided before aliquoting the urine. If multiple 24hr containers for the same test was provided – combine both 24hr containers to measure total volumes before aliquoting	All 24hr urine tests	Multiple	URB	Clear bottle (60ml) or Amber bottle (60ml) 24 hr. urine tests requiring protection from light – aliquot into amber bottle 
Urine Culture Specimens for urine culture must be submitted in Grey top tube – all other containers will be canceled	32	Multiple	MIC	Grey top – Boric Acid Tube 
Urine Immunology Monogam Initial, Random Monogam Initial, Timed Monogam Follow Up, Random Monogam Follow Up, Timed Protein Electrophoresis	UMIR UMIH UMFR UMFH UER	4634 4737 4738 4739 4518	URN/URB URN/URB URN/URB URN/URB URM/URB	Clear bottle (60ml) and Clear Top Discard Tube 

Department of Pathology and Laboratory Medicine

Urine Collection Device Grid



Version 08/2025

Test Name	SoftLab Code	Epic Code	Label Code	Collection Device
Urine Chemistry				Clear Discard Tube – No Preservative
Amino Acids, Quantitative, Random	UAARZ	5306	URN	
Amylase, Random	UAMYR	9551	URN	
Amylase/Creatinine Ratio	AMYCR	5568	URN	
Calcium, Random	UCAR	9510	URN	
Calcium/Creatine Ratio	UCCR	9525	URN	
Creatinine, Random	UCRR	5114	URN	
Electrolytes	ULYTH	9550	URN	
Microalbumin, Random	UMALR	9472	URN	
Microalbumin/Creatine Ratio	UMCR	4580	URN	
Myoglobin	MYOGX	5552	URN	
Protein, Total, Random	UTPR	5609	URN	
Protein/Creatinine Ratio	UPCR	5087	URN	
Urea Nitrogen, Random	UUNR	9506	URN	
Uric Acid, Random	UUAR	9516	URN	
Urine Pregnancy	UPREG	91	URN	Clear Discard Tube – No Preservative
				
Molecular				Alinity Multi-Collect Specimen Kit
Chlamydia DNA by PCR	CHLPC	1026	MST	
Chlamydia/GC DNA by PCR	CHGC	1114	MST	

Department of Pathology and Laboratory Medicine

Urine Collection Device Grid

Version 08/2025

Test Name	SoftLab Code	Epic Code	Label Code	Collection Device
Cytology	CYU			Cytolyte Container (30ml) or Non-Sterile Urine Cup 
Drug Screen				
Comprehensive Drug Analysis	UPDSZ	2535	URN/URB	Clear bottle (60ml) and Clear Top Discard Tube 
Drug Abuse Screen, 6 panel w/Confirmation	UDS6	5908	URN/URB	
Drug Abuse Screen, 6 panel w/o Confirmation	UDS6	9602	URN/URB	
Drug Abuse Screen, 8 panel w/Confirmation	UDS8	5911	URN/URB	
Drug Abuse Screen, 8 panel w/o Confirmation	UDS8	9601	URN/URB	