

<b>TEST</b>	<b>CONTAINER TYPE</b>
Giardia Antigen Stool by EIA Giardia/Cryptosporidium Antigen Stool by EIA Stool Culture (Routine), includes E.coli 0157:H7 Vibrio Stool Culture Yersinia Stool Culture	<p style="text-align: center;"><b>C&amp;S Fecal Swab</b></p> 
Enteric Pathogens Panel by NAAT Giardia Antigen Stool by EIA Giardia/Cryptosporidium Antigen Stool by EIA Stool Culture (Routine), includes E.coli 0157:H7 Vibrio Stool Culture Yersinia Stool Culture	<p style="text-align: center;"><b>C&amp;S PARA PAK</b></p> 
Giardia Antigen Stool by EIA Giardia/Cryptosporidium Antigen Stool by EIA Ova and Parasites Examination	<p style="text-align: center;"><b>PVA and Formalin Vials</b></p> 
Clostridium difficile Antigen/Toxin Helicobacter Pylori Antigen Stool Occult Blood - Stool Parasite ID (Identification) Rotavirus Stool pH Vibrio Stool Culture WBCs Stool Yersinia Stool Culture	<p style="text-align: center;"><b>Stool Cup</b></p> 

<p>Chlamydia Trachomatis by NAAT  N. Gonorrhoeae by NAAT  Trichomonas Vaginalis NAAT</p>	<p><b>Xpert CT/NG Vaginal Endocervical Specimen Collection Kit</b></p> 
<p>Fluid Culture w/Anaerobes and Smear (Body Fluid Culture)  Wound Culture w/Anaerobes and Smear (Any Wound Aspirate)</p>	<p><b>Sterile Syringe</b></p> 
<p>CSF Culture w/Smear  Cytomegalovirus (CMV) Culture  Herpes Simplex Virus (HSV) 1&amp;2 DNA by NAAT (CSF)  Meningitis/Encephalitis Pathogen Panel by NAAT CSF</p>	<p><b>CSF Vials</b></p> 
<p>2019 Coronavirus SARS-CoV-2 FL  Herpes Simplex Virus (HSV) 1&amp;2 DNA by NAAT (NOT CSF)  Influenza A,B by NAAT  Respiratory Syncytial Virus (RSV) by NAAT  Respiratory Pathogens Panel by NAAT</p>	<p><b>Viral-Universal Transport Media</b></p> 
<p>Bronchial Culture  Chlamydia Trachomatis by NAAT Urine  Catheter Tip Culture  Cystic Fibrosis Respiratory Culture w/Smear (CF Patients Only)  Cytomegalovirus (CMV) Culture  N. Gonorrhoeae by NAAT Urine  Respiratory Culture w/Smear - Sputum  Tissue Culture w/Anaerobes &amp; Smear  Trichomonas Vaginalis by NAAT (Urine)</p>	<p><b>Sterile Container</b></p> 

Bordetella Pertussis by NAAT  
Cystic Fibrosis Respiratory Culture w/Smear (CF Patients Only)  
Ear, Eye, or Skin Culture w/Smear  
GC Culture w/o Smear  
Genital Culture w/Smear  
Group A Streptococcus by NAAT (Throat)  
Group B Streptococcus by NAAT  
MRSA Screen by NAAT  
Mycoplasma Pneumonia by NAAT  
Nasal MSSA / MRSA by NAAT  
Strep Screen Confirmation (Beta Hemolytic Streptococcus Culture (Throat) (Group A)  
Throat or Nasopharyngeal Culture  
VRE by NAAT (Vancomycin Resistant Enterococcus)  
Wound Culture w/Anaerobes and Smear

**E-Swab**



Cytomegalovirus (CMV) Qualitative NAAT (Saliva)

**Copan 553C collection device**

