

# Guide for Microbiology Collection

## GC & Chlamydia by PCR

\*Female Only  
Endocervical or  
patient collected vaginal

Collect & place swab from  
Xpert CT/NG  
Vaginal Endocervical  
Specimen Collection Kit  
into **pink top** transport  
reagent tube



## GC & Chlamydia Urine by PCR

\*Female or Male  
Random Urine

Collect urine in sterile, leak  
proof container



## Viruses by PCR

HSV  
VZV  
CMV  
Influenza A&B  
RSV by PCR  
RPP by NAAT

Collect & place swab from  
UTM (Universal Transport  
Media) Kit into **red top**  
transport reagent tube



## Bacterial Culture

Copan Liquid Amies  
Elution (eSwab)  
Collection & Transport  
System Kit

Collect & place swab  
from eSwab Kit &  
break off swab into  
transport tube

